



Требования к биоматериалу

Краткая таблица для подготовки и передачи

Материал	Основные требования
Свежая ткань	Перед забором согласуйте среду, температуру и время доставки. Не замораживать и не промывать без инструкции.
Фиксированная ткань	Герметичный контейнер, 10% нейтральный забуференный формалин. Для ER/PR/HER2 молочной железы обычно 6-72 ч фиксации.
Парафиновые блоки	Жёсткая упаковка, комнатная температура, защита от жары и солнца.
Готовые стекла	Маркированные стекла в жёстком боксе; приложить заключение и опись.
Цитологический мазок	Правильная фиксация по типу исследования; защитить от царапин и поломки.
Жидкостная цитология	Только согласованный контейнер и консервант.
Кровь / костный мозг	Пробирка, объём, температура и сроки зависят от теста; согласовать до забора.
Декальцинированная ткань	Указать метод декальцинации: сильная кислота может ухудшить ИГХ и молекулярные тесты.

Важно: точные условия зависят от конкретного теста. До забора нового материала подтвердите требования у CellDig.



Biomaterial talablari

Tayyorlash va topshirish bo'yicha qisqa jadval

Material	Asosiy talablar
Yangi to'qima	Olishdan oldin muhit, harorat va yetkazish vaqtini kelishing. Ko'rsatmasiz muzlatmang yoki yuvmang.
Fiksatsiyalangan to'qima	Germetik idish, 10% neytral buferlangan formalin. Ko'krak ER/PR/HER2 uchun odatda 6-72 soat.
Parafin bloklar	Qattiq qadoq, xona harorati, issiqlik va quyoshdan himoya.
Tayyor slaydlar	Belgilangan slaydlar qattiq boksda; xulosa va ro'yxatni qo'shing.
Sitologik surtma	Tekshiruv turiga mos fiksatsiya; tiralish va sinishdan himoya.
Suyuqlik sitologiyasi	Faqat kelishilgan idish va konservant.
Qon / suyak ko'migi	Probirka, hajm, harorat va muddat testga bog'liq; olishdan oldin kelishing.
Dekalsinatsiyalangan to'qima	Usulni ko'rsating: kuchli kislota IHK va molekulyar testlarga xalaqit berishi mumkin.

Muhim: aniq shartlar testga bog'liq. Yangi material olishdan oldin CellDig bilan tasdiqlang.



Specimen requirements

Quick preparation and transport table

Material	Key requirements
Fresh tissue	Confirm medium, temperature and transport time before collection. Do not freeze or rinse unless instructed.
Fixed tissue	Leakproof container with 10% neutral buffered formalin. Breast ER/PR/HER2 specimens are generally fixed for 6-72 hours.
Paraffin blocks	Rigid packaging at room temperature; protect from heat and sunlight.
Prepared slides	Labelled slides in a rigid box; include report and inventory.
Cytology smear	Use fixation appropriate to the test; protect from scratches and breakage.
Liquid-based cytology	Use only the agreed container and preservative.
Blood / bone marrow	Tube, volume, temperature and timing are test-specific; confirm before collection.
Decalcified tissue	State the decalcification method; strong acid may impair IHC and molecular testing.

Important: exact conditions depend on the test. Confirm with CellDig before collecting a new specimen.